**geemarc®** Wake Up Call  
Hear the Difference

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



# DESCRIPTION

En

(\*) Classified TVN-3 according to EN60950 standard.

## 1- CONNECTING THE TELEPHONE

- 1 - Connect the PTT line cord.
- 2 - At your option you may connect your telephone in the doubler plug.
- 3 - Connect the PTT doubler plug in the wall socket.
- 4 - Connect the adaptor line cord.
- 5 - Connect the adaptor in the wall socket. (\*\*)

(\*\*)Classified "hazardous voltage" according to EN60950 standard.

**Electrical connection:** The apparatus is designed to operate at DC9V, 650mA. The power adaptor supplied operate from AC100 – 240V, 50/60Hz input.

The apparatus does not incorporate an integral power on/off switch. To disconnect the power, either switch off supply at the mains power socket or unplug the DC adaptor. When installing the apparatus, ensure that the mains power socket is readily accessible.

**Telephone connection:** Voltages present on the telecommunication network are classified TNV-3 (Telecommunication Network Voltage) according to the EN60950 standard.

The lightning flash with arrowhead symbol, within and equilateral triangle, is intended to alert the user to the presence of uninsulated "dangerous voltage" within the product's enclosure that may be of sufficient magnitude to constitute a risk of electric shock to persons.

**CAUTION:**

TO REDUCE THE RISK OF ELECTRICAL SHOCK, DO NOT REMOVE THE COVER (OR BACK) FROM THE PRODUCT. NO USER SERVICEABLE PARTS INSIDE. REFER SERVICING TO YOUR SONIC ALERT DEALER.

The exclamation mark within an equilateral triangle is intended to alert the user to the presence of important operating and maintenance (servicing) instructions in the literature accompanying the appliance.

**FCC Part 15 statement:**

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: 

- (1) this device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

**This product meets the applicable Industry Canada technical specifications.**

**The Ringer Equivalence Number is an indication of the maximum number of terminals allowed to be connected to a telephone interface. The termination on an interface may consist of any combination of devices subject only to the requirement that the sum of the Ringer Equivalence Numbers of all the devices does not exceed five.**

## 2- SETTING THE TIME ( )

---

- Press and hold  (set clock) button.
- Press **H** (hour) and **M** (minute) until correct hour and minutes appear on the display.
- Release  (set clock) button.
- For p.m. hours, the p.m. indicator will appear on the display near the inscription **PM**.

### 12 / 24 hour format

User can select which hour format to use. To select:

- Press and hold the **H** (hour) button for 6 seconds to enter the hour selection mode, the display shows “12H” or “24H”, then release the hour button.
- Press the **H** (hour) button to toggle between 12 or 24 hours.

**Note:** If no selection made within 4 seconds, the last showed format will be chosen, and the clock will return to normal time mode.

## 3- SETTING THE ALARM TIME ( )

---

- Press and hold  (set alarm) button.
- Press **H** (hour) and **M** (minute) until correct hour and minutes appear on the display.
- Release  (set alarm) button.
- Press , **H** and **M** to reset the alarm clock to 12.

## 4- SETTING THE ALARM MODE

---

Put the function switch on **A+T** position.

Choose which alarm mode you wish to be awakened to by moving the alarm mode switch to one of the following selections:


**V** : only shaker activated

**V+S** : both sound and shaker alarm activated

**V+F** : both strobe light and shaker alarm activated

**V+F+S**: strobe light, sound and shaker alarm activated

**Note:** if user doesn't use the **V** (shaker function), please unplug the shaker plug from the clock.

When the alarm is set to one of the alarm modes, the alarm indicator will appear on the display near the inscription **AL**.

The table below shows how the alarm modes indicate an incoming call from the alarm clock.

ALARM MODE	ACTION	ALARM CLOCK	INCOMING CALL
SOUND	Buzzer	With pauses	With pauses
SHAKER	Vibrates	With pauses	With pauses
STROBE LIGHT	Flashes	With pauses	With pauses

The default alarm duration time is 30 minutes. User can adjust the alarm duration time from **1-59** minutes.

To adjust the alarm duration:

- Press and hold the set alarm button  for 6 seconds to enter adjust mode. The clock show "30", release the set alarm button.
- Press the **M** button to adjust the alarm duration time to the desired setting from 1 minute to 59 minutes.
- This new setting will become the new default, until changed by user.
- If no action is take within 3 seconds, then the clock will return to normal time mode.

## 5- SNOOZE KEY ( )

---

When the alarm mode you have chosen is switched on the alarm time setting, pressing the  button will stop the alarm mode for approximately 9 minutes. The alarm will then switch the alarm mode on again for a further 9 minutes. This sequence can be carried out for up to 59 minutes after which the alarm mode will be automatically switched off.

User can adjust the snooze repeat time from **1-30** minutes.

- Press and hold the **SNOOZE** button for 4 seconds to enter the snooze time adjust mode. When the clock show “9” release the Snooze button.
- Press the **M** button to adjust the snooze time to the desired setting from 1 minute to 30 minutes.

This new setting will become the new default, until changed by user. If no action is taken within 3 seconds, then the clock will return to normal time mode.

## 6- SWITCHING THE ALARM OFF FOR 24H. (OFF)

---

Press the alarm off button (**OFF**), the alarm will then operate at the same time the next day. The alarm if not switched off will continue to sound for up to 59 minutes and then be automatically switched off.

## 7- RINGER On/Off SWITCH ( A+T / OFF / T )

---

- **A+T** (alarm + telephone): for the alarm and telephone ring position.
- **OFF**: alarm and ring off.
- **T** (telephone): for the telephone ring position.

## 8- TELEPHONE OPERATION

---

- Ensure that the telephone has been set up in accordance with the telephone userguide and has been connected as shown in the installation section.
- Select the ringer On/Off switch on the **A+T** or **T** position for the ring.
- Choose which of the alarm modes you wish to alert you of your incoming telephone calls (see page 4: **V**, **V+S**, **V+F**, **V+S+F**).
- When the alarm mode selected indicates there is an incoming telephone call, lift the telephone handset and answer the call. When you have finished the call replace the handset and the alarm will be ready to indicate the next call.

## 9- Dimmer Function ( )

---

There is a 5-level brightness control for the LED display. In normal time model, press the **Snooze** button once to change the display brightness. The brightness will cycle from bright to dim then from dim to bright (the cycle repeats)

**NOTE** : The dimmer function will not work when the clock is in snooze mode.

## **10- VOLUME ( VOLUME )**

---

Sets the volume of the ring and the alarm sound.

## **11- TONE ( TONE )**

---

Sets the tone of the alarm sound.

## **12- SHAKER SWITCH ( + ((O)) - )**

---

Sets the motion of the shaker to Lo (-) or Hi (+).

## **13- BATTERY BACK- UP**

---

It's not necessary to insert a battery for the clock to function. With battery back-up installed, the clock will keep the time and alarm settings in the event of a power cut or the DC adaptor being unplugged.

Check the battery condition after one year, replacing old battery. A leaking battery can destroy your clock.

- This equipment complies with Part 68 of the FCC rules and the requirements adopted by the ACTA. On the bottom of this equipment is a label that contains, among other information, a product identifier in the format US: AAAEQ##TXXXX. If requested, this number must be provided to the telephone company.

- All applicable certification jacks Universal Service Order Codes (USOC) for the equipment is provided (i.e. RJ11C) in the package with each piece of approved terminal equipment.

- A plug and jack used to connect this equipment to the premises wiring and telephone network must comply with the applicable FCC Part 68 rules and requirement adopted by the ACTA. A compliant telephone cord and modular plug is provided with this product. It is designed to be connected to a compatible modular jack that is also compliant. See installation instructions for details.

- The REN is useful to determine the quantity of devices you may connect to your telephone line and still have all of those devices ring when your telephone number is called. In most, but not all areas, the sum of the REN's of all devices connected to one line should not exceed five (5.0). To be certain of the number of devices you may connect to your line, as determined by the REN, you should contact your local telephone company to determine the maximum REN for your calling area.

- If this equipment causes harm to the telephone network, the telephone company will notify you in advance that temporary discontinuance of service may be required.

But if advance notice isn't practical, the telephone company will notify the customer as soon as possible.

Also, you will be advised of your rights to file a complaint with the FCC if you believe it is necessary.

- The telephone company may make changes in its facilities, equipment, operation or procedures that could affect the operation of the equipment. If this happens the telephone company will provide advance notice in order for you to make necessary modifications to maintain uninterrupted service.

- This equipment may not be used on coin service provided by the telephone company. Connection to party lines is subject to state tariffs. This equipment is hearing aid compatible.

**Warning:** Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

**NOTE:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
  - Increase the separation between the equipment and receiver.
  - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
  - Consult the dealer or an experienced radio TV technician for help.
- If trouble is experienced with this equipment, for repair or warranty information, please contact our customer service staff at: **1-888-864-2446 (RMA@sonicalert.info)**

This class B device applicable Industry Canada technical specifications.

If the equipment is causing harm to the telephone network, the telephone company may request that you disconnect the equipment until the problem is resolved.

- Please follow instructions for repairing if any; otherwise do not alter or repair any part of device except as specified.

**Opening the equipment or any attempt to perform repairs will void the warranty. For service or repairs, call 1-888-864-2446. (RMA@sonicalert.info)**

- Connection to party line service is subject to state tariffs. Contact the state public utility commission, public service commission or corporation for information.
- If your home has specially wired alarm equipment connected to the telephone line, ensure the installation of this telephone equipment does not disable your alarm equipment. If you have questions about what will disable alarm equipment, consult your telephone company or a qualified installer.

Sonic Alert warrants the product against any defect in materials or workmanship for the period of one year from the date of purchase.

Should you experience a problem, contact our customer service department. Be sure to save your sales receipt as proof of purchase date should you need warranty service.

Within a period of one year from purchase date, Sonic Alert will repair or replace (our discretion) your product at no cost, if a defect in materials or workmanship is found.

If we elect to replace your product, we may replace it with a new or reconditioned product of the same or similar design.

Repair or replacement will be warranted for a period of 90 days or the original time on the original warranty, whichever is longer.

The warranty does not cover accidents, negligence or breakage to any parts. This includes shipping damage, failure to follow instructions, misuse, fire, floods, use of incompatible accessories, Acts of God or failure in your product service carrier's line service. The product must not be tampered with or taken apart by anyone who is not an authorized Sonic Alert representative. Tampering with the product will void any written or implied warranties. Sonic Alert shall not be responsible for loss of time, inconvenience, property damage caused by your product or any other accidental or consequential damages.

Warranty service is available only with proof of purchase.

Simply send the Wake Up Call(postpaid) and a copy of your sales slip as proof of purchase to:

**Sonic Alert Warranty Center 1081 West Innovation  
Dr. Kearney, MO 64060** (RMA@sonicalert.info)

**IMPORTANT : YOUR RECEIPT IS PART OF YOUR  
WARRANTY AND MUST BE RETAINED AND  
PRODUCED IN THE EVENT OF A WARRANTY CLAIM.**

# SonicAlert

Sonic Alert Warranty Center  
1081 West Innovation Dr.  
Kearney, MO 64060

Sales: 248-577-5400 Fax: 248-577-5433  
Customer Service & Warranty: 1-888-864-2446  
RMA@sonicalert.info

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